

LED Turbo 2100mw/cm²

3 Second, Cordless LED - Super Fast Curing Light

Features

- Super High Output (2,100 mw/cm²) for super fast curing (3 seconds up to 3 mm)
- 3 curing modes - *Standard, Ramp and Boost.*
- High capacity battery - best battery in the business (2 year warranty) get over 1 hour of curing time on full battery charge (that is 1,200 three second cures on one battery charge)
- Fast charging base with built-in testing device. Get a 90% battery charge on 1 hour of charging time.
- Ergonomic gun type design
- High quality glass light guides, 8mm turbo tip included. Optional 11mm, 8mm standard or 3mm tacking tip available.



Optional Light Guides



Complete Kit Includes:

Curing Light, 8mm Light Guide, Light Shield, Charging Base, Power Transformer.

1 year warranty, 2 year warranty on battery

Item # VLEDT Retail \$ 700.00 ea

Special
\$595⁰⁰

Offer Expires 12/31/2010

Vector R & D Inc.

6824 19th St. # 230

University Place, WA 98466

877-883-7455 www.vectorusa.net

Available from your local dental supply dealer

