

# Lil' BEAVER 2.0

## MAGNETOSTRICTIVE Ultrasonic Scaling Unit



- Accepts both 25K (long) and 30K (short) ultrasonic inserts - no adapter or switch, unit automatically adjusts to accommodate any style, brand or frequency.
- Convenient front-mounted water control with large easy-to-adjust dial - easy even with gloves.
- Turbo boost, 2 position foot pedal - press lightly for regular dial power, press firmly for TURBO! *Heavy calculus watch out!*
- Digital power adjustment with LED lit power display panel.
- Front mounted main power switch - shutting the unit off when left unattended just got a whole lot easier.
- All aluminum housing for improved durability. No more plastic to stain, crack or break
- Lies flat or mounts vertically in provided vertical mounting stand for optimal counter space saving. Fits almost anywhere!

**Includes vertical  
mount base!**



### Complete Kit Includes:

Little Beaver Scaling Unit, Power Cord, 2 Position Turbo Boost Foot Pedal, 1/4" Water Line w/Quick-Disconnect, Vertical Mount Stand, Handpiece Holder Mount, Autoclavable Handpiece Sheath. 1 year warranty.

SC-001 ..... Retail \$ 500.00 ea

**Vector R & D Inc.**  
6824 19th St. # 230  
University Place, WA 98466  
877-883-7455 [www.vectorusa.net](http://www.vectorusa.net)

Available from your local dental supply dealer

