

Little Beaver

Magnetostrictive Ultrasonic Scaling Unit

By Vector R & D Inc.

Dual Frequency Auto-adjust 25K/30K Ultrasonic Scaler

The beaver is well known as one of the hardest working and most clever creatures in the animal kingdom. Just like its name sake, the Little Beaver Scaler is designed for hard working durability and packed with clever features like:

- 25K/30K auto-frequency adjustment that allows you to use all your 25K and/or 30K "cavitron™" type inserts.
- Two position "Turbo Boost" foot pedal. Press lightly to begin scaling. Press more firmly and the scaler kicks up to a higher power to attack an area of heavier calculus. All this is accomplished with the foot control, so there is no need to stop scaling and adjust the power dial.
- Removable, sterilizable handpiece sheath. The outer housing of the handpiece is removable and can be sterilized. Eliminates the need for cumbersome bags and plastic wraps.



Compare and save and you will choose a "Beaver" over a "Bobcat" any day!

Cool and clever features!



Easy and precise water control adjustment. True needle valve water control provides precise and reliable water adjustment.



Removable, autoclavable handpiece sheath. Eliminates the need for cumbersome bags and barrier wraps. Our competitors require that you sterilize the entire handpiece, while the Little Beaver only requires sterilization of the the outer housing keeping all the electrical components out of the autoclave, thus eliminating repairs and/or costly handpiece replacements.

Automatic 25K 30K frequency adjustment!
The Little Beaver scaler automatically adjusts the scaling frequency based on the type of insert placed in the handpiece. Use any and all of your existing inserts! Any brand, any manufacturer, any frequency, the Little Beaver takes them all!

Little Beaver Ultrasonic Scaler with Auto-adjust 25K/30K tip compatability by Vector R & D Inc.

Kit Includes: Unit, Foot Control, Power Cord, Handpiece, 1 Autoclavable Removable Handpiece Sheath, Water Tube with Quick Disconnect.

1 year warranty. Extra sterilizable handpiece sheaths sold separately.

Scaler complete kit SC-001 Sug. Retail \$ 500.00 ea

Extra handpiece sheaths SC-050 Sug. Retail \$ 35.00 ea

Complete Kit:



Note: scaling tip not included

Vector R & D Inc. 6824 19th St. #230 University Place, WA 98466
Telephone 877-883-7455, Fax 253-564-4927 web www.vectorusa.net
Available from your dental distributor