

Vector R&D Presents:

THE

Turbo-torque

America's Best Value Low Speed Motor System!

Compare!
Get more for your money
with Vector Turbo Torque!

Push Button
Head Included!



Versatile E-Type Connection.

E-Type low speeds are quickly becoming the "standard" for the industry. Easily interchanges with most electric, surgery and endo systems. With our competitors you are stuck with the attachments that they build to fit their handpiece. With Vector you can choose from not only our high quality attachments but from any manufacture who manufactures E-type attachments.

Top Quality, High Torque Motor

While our attachments are some of the highest quality in the dental industry, the real corner stone of this set is the motor. This smooth and quiet motor is designed to last for years and years. A real proven performer with consistent high torque throughout the entire speed range. *Do not let the price fool you this is a top quality motor!*



20,000 RPM our most popular motor Set

Our most popular 20,000 rpm design is perfectly suited for operative work and acrylic adjustment as well as high speed polishing of the newer composites. Add a 4:1 straight nose cone (next page) and the handpiece is transformed to the perfect "prophy" motor. Comparable packages for our competitors will cost well over \$ 1,000.00. The Turbo Torque is truly value without compromise!



VECTOR
RESEARCH & DEVELOPMENT INC.

6824 19th St. #230, University Place, WA 98466
Tel. 253-564-5084 - Fax 253-564-4927

Toll Free 1-877-88DRILL - WWW.VECTORUSA.NET

Dealer
inquires
welcome!

Vector Turbo Torque Product Specifications

Vector Turbo Torque

Compare to other motor systems and **SAVE!**

Low Speed Deluxe Motor Set

Complete 20,000 RPM motor set. Includes 20,000 rpm motor, 1:1 contra angle attachment, push button latch type contra angle head, and a 1:1 straight nose cone. Add additional parts and options from the selection below to complete the full functional range of this versatility motor system. Connects to standard 4 hole tubing. Motor will accept any standard ISO E-Type contra angle. Contra angle will accept any NSK fitment contra angle head.

20,000 Rpm Vector # (TT-20K)

Complete motor set!



Turbo Torque Attachments

1:1 No Reduction Straight Nose Cone



4:1 Gear Reduction Straight Nose Cone



High-Performance Attachments!

1:1 Contra Angle



4:1 Contra Angle



16:1 Contra Angle



FG Head

Latch Head

TURBO TORQUE

Straight Nose Cone Attachment 1:1

Fits all standard E-Type connections. Will accept any 3/16" shaft lab bur or "straight handpiece bur or any "doriot" style prophylaxis or contra angle.

1:1 Straight attachment Vector # (TT-SN1)

TURBO TORQUE

Straight Nose Cone Attachment 4:1

4:1 Reduction gear generates 5,000 rpm from a 20,000 rpm motor. Fits all standard E-Type connections. Will accept any 3/16" shaft lab bur or "straight handpiece" bur or any "doriot" style prophylaxis or contra angle.

4:1 Reduction straight attachment Vector # (TT-SN4)

TURBO TORQUE

1:1 Contra Angle Attachment

Standard no reduction contra angle. Accepts all NSK fitment contra angle heads and attachments. Will fit any standard E-type motor.

1:1 Contra angle attachment Vector # (TT-EC)

TURBO TORQUE

4:1 Reduction Contra Angle Attachment

Fits all standard E-Type connections. Will accept any "NSK" threaded contra angle head.

4:1 Contra angle attachment Vector # (TT-E4R)

16:1 Reduction Contra Angle Attachment

Fits all standard E-Type connections. Will accept any "NSK" threaded contra angle head.

16:1 Contra angle attachment Vector # (TT-E16R)

Push Button Head Attachments

For any "NSK" threaded contra angle attachment. Will not fit "Midwest" contra angles.

Friction grip head attachment

Vector # (TT-LAP)

RA latch type head attachment

Vector # (TT-FGP)

Don't Forget the Vector Spray!



Dealer inquires welcome!



VECTOR
RESEARCH & DEVELOPMENT INC.

6824 19th St. #230, University Place, WA 98466
Tel. 253-564-5084 - Fax 253-564-4927

Toll Free 1-877-88DRILL - WWW.VECTORUSA.NET