



Integrated plug-in monitor, using the A series of plug-in modules, plug and play

8 inch color LCD screen, multi-channel waveform display

The module is complete, can be any combination: support EMS module box (optional 3/5/12 guide ECG, RESP, NIBP, TEMP, IBP, SpO2), CO2, AG, ICG and other modular box

Flexible input function: buttons, knobs, touch screen, mouse, keyboard and so on

Free combination of multi-screen, multi-interface display
Configurable memory card, trend data is saved, Configurable recorder, record monitoring data and waveforms



3.5 inch HD color LCD display, multi-channel waveform display

Full-screen touch, quick button operation, to facilitate the operation of clinical staff

Lightweight and portable, with onboard power interface.

With gravity sensing, so that medical staff at any angle to observe the patient's monitoring data and waveform

Multi-parameter monitoring module EMS is the carrier of the monitoring data, compatible with all the A series of monitors, to achieve continuous monitoring.



MT-200SET	Complete Entitled C02 8" Monitor set (MT-200) + EMS 3 module + ETCO2 Module (MT-220)	\$ 5,450.00
MT-200	8" TFT LCD wide Screen, EMS3 Module 3.5" (Standard config + BLT SpO2 + 3/5 lead ECG,) resolution 800x600, Arrhythmia/ ST segment analysis, 168 hours graphic / list trending, VGA display interface, 2 USB ports (peripheral devices interface)	\$ 3,695.00
MT-100	EMS3 Module 3.5" (Standard config + BLT SpO2 + 3/5 lead ECG), with Docking Station	\$ 1,990.00
MT-110	EMS 3 docking / charging station	\$ 450.00
MT-220	Respironics LoFlo EtCO2 sensor with module box	\$ 2,150.00
MT-200M	8" Monitor Only MT-200 Monitor Only Upgrade kit for MT-100 to make 8" monitor	\$ 1,800.00