

45° Highspeed - No Air Jet!



New! Vector ECO 45

The new Vector ECO series 45 is an excellent value and great combination of durability and performance. No chip-air into the surgical sight reduces the risk of air embolism when compared to a standard air driven highspeed handpiece. 45 degree angle design improves access to the 3rd molar region. Quick-release swivel connector easily removes for sterilization between patients. Medium head provides enough torque and power without obstructing the view of the surgical sight when compared to large head air-driven handpieces. Vector engineering and design provides quiet and efficient operation with plenty of power to get the job done.



- Excellent torque and cutting efficiency!
- Improved bur holding power!
- Direct replacement for NSK CH45 style highsPEEDS!
- Single water spray - no chip air!



Vector ECO 45

High-performance air-driven 45 degree, no-chip-air, highspeed handpiece. With quick connect swivel connection included, connects to 4 hole non-optic tubings. 6 month warranty on turbine, 1 year on handpiece. Do not autoclave the swivel connector, handpiece only.

Includes: Handpiece, 4 hole swivel connector, lubrication adapter, cap wrench. 6 month warranty on turbine, 1 year on handpiece.

E-45-4 Retail \$ 400.00

Accepts NSK # NCH-SU03 or Vector # TNCS

TNCS Retail \$ 175.00