


2810 69th Ave West - University Place, WA 98466
www.vectorusa.net 877-88DRILL (877-883-7455), Fax 253-564-4927

All New MD Technologies Implant Motor Set!

Complete set includes: Control box, Full-function- variable speed foot pedal, 40,000 rpm motor with led light, 20:1 reduction fiber-optic contra angle, integrated water pump with sterile water tube. Everything you need to get started! Just get your implant starter kit and go! Fully adjustable and programmable works with every brand make or model of implant.



Programable Full-function control box



FEATURES:

- a) Adjustable pre-programmed gear ratios
- b) Adjustable water flow control
- c) forward / reverse adjustment
- d) Light on / off
- e) Program setting / save mode
- f) Speed , torque and program selection.
- g) on-board peristaltic pump for sterile water delivery.

MD-Tech motor complete set: Control box, foot control, 20:1 handpiece, handpiece stand, tubing set, accessories.

- MT-IO20.....Reg. \$2,700⁰⁰
- 20:1 Optic Implant Contra Angle
- MT-IO20.....Reg. \$695⁰⁰
- Implant Tubing Set
- MT-IO20.....Reg. \$20⁰⁰