

LET'S DO IMPLANTS!

Complete implant motor set! All you need is the implants!



Complete Set:
Control Box + Foot Control + Motor + 20:1 Angle!

Now with LED
Lighted Motor and
Fiberoptic Angle

Only
\$2,495⁰⁰



20:1 Optic
Handpiece!

Full function
variable speed
foot pedal!

All New MD Technologies Implant Motor Set!

Complete set includes: Control box, Full-function- variable speed foot pedal, 40,000 rpm motor with led light, 20:1 reduction fiber-optic contra angle, integrated water pump with sterile water tube.

- Everything you need to get started!
- Get your implant starter kit from your favorite implant manufacturer and go!
- Fully adjustable and programmable.
- Works with every brand make or model of implant.

CONTROL BOX FEATURES:

Lighted LCD Display

Surgical & Implant settings!

- 1) Adjustable pre-programmed gear ratios
20:1, 27:1, 32:1, 1:1, 1:2, 16:1
- 2) Adjustable irrigant water flow control
- 3) Forward / reverse adjustment
- 4) Light on / off
- 5) Program setting / save mode



Adjustable

- Speed
- Torque
- Program selection
- Water Flow
- LED Light on/off

20:1 Angle
Speed Range 2,000 - 10rpm
70nm - 5nm

Complete Implant System!



MD-Tech Complete Motor Set:

Includes: Control box, foot control, 20:1 Optic handpiece, multi-function variable speed foot control, handpiece stand, tubing set, accessories.

MT-ISE..... Reg. \$3,200⁰⁰ Special **\$2,495⁰⁰**

20:1 Optic Implant Contra Angle
MT-IO20A..... Reg. \$950⁰⁰ Special **\$695⁰⁰**

Implant Tubing Set
MT-ISE-TB..... Reg. \$100⁰⁰ Special **\$70⁰⁰**

30 Contact Your Dental Dealer To Place an Order